Safety instructions and general instructions for BEGO investment materials

BEGO is an official co-partner of the German Olympic teams for **Salt Lake City in 2002** and **Athens in 2004**

Safety instructions

- Investment materials contain quartz. **Do not inhale dust!** Risk of damage to lungs (silicosis, lung cancer).
Recommendation: Use gas mask type FFP 2 – EN 149:2001. Cut open bag with scissors and avoid generation of dust when filling into mixing bowl. Rinse out empty bag with water before crumpling up. Remove dust at workplace only when **moist**. To avoid dust during deflasking, place the mould in water after it has completely cooled down after casting until it is thoroughly moistened. Use suction extraction system with fine dust filter when blasting.
- Remove dried-up mixing liquid only when **moist** (fine dust). Seal bottle securely after use.
- Suction-extract escaping furnace gases during preheating and discharge outdoors!
- Only with shock heat investment materials: **Risk of injury during shock heating!** Place all moulds in the furnace within 10 seconds – then keep the furnace door closed for 15 minutes!

General instructions

- The ideal processing temperature for BEGO investment materials is **20 °C**. To keep it constant at a higher ambient temperature, use a temperature control cabinet if necessary. At a high ambient temperature also place the mixing bowl and mixing bowl in the temperature control cabinet.
- Observe shelf life dates of investment material and mixing liquid! Do not use without checking after the specified date. Do not bring phosphate-bonded investment material into contact with plaster or investment material containing plaster.
- Do not use crystallized mixing liquid any longer. The higher the concentration of mixing liquid, the greater the expansion. **A practical tip: Keep spare bottles!** Fill BegoSol® bottle up to desired % mark with BegoSol® and top up to 100 % with distilled water. Label this spare bottle with % data. Advantage: The mixing liquid is available with the desired concentration and can be measured in one operation.
- Data on the concentration of the mixing liquid are based on our own experience and tests and are regarded as reference values. Work results may be influenced by such parameters as temperature, mixing intensity, investment system, etc.

Application-related recommendations made by us, whether given verbally, in writing or by way of practical instructions, are based on our own experience and tests and may therefore be regarded only as general guidelines.

Our products are subject to continuous development. We therefore reserve the right to change the design and composition of our products.

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Precision casting investment material for partial dentures
– shock or conventional heating

Safety instructions


Please read and follow the instructions in the insert
“Safety instructions and general instructions
for BEGO investment materials”!

BEGO is an official co-partner of
the German Olympic teams for
Athens in 2004

WiroFine can be heated up rapidly (“shock heat”) or conventionally.


The possible insertion temperatures are: approx. 20 °C / 70 °F (conventional) or 700 °C to 1000 °C / 1290 °F to 1830 °F (shock heat).

General instructions

- Liquid for shock and conventional preheating:
BegoSol® K (**Frost-sensitive!** Storage and transport temperature: + 5 °C to + 35 °C / 41 °F to 9
- Before mixing, rinse out the clean mixing bowl with water and wipe off.
Mixing bowls that are not clean or are dry withdraw moisture from the investment material!
- Processing width 20 °C / 70 °F: approx. 3.5 minutes
23 °C / 73 °F: approx. 3.0 minutes
- First put in liquid and add powder, mix thoroughly with a spatula for at least **15 seconds**.
Then mix in a mixing unit, such as *EasyMix*, for 60 seconds under vacuum conditions, as far as possible.
(Processing without mixing unit: mix for **2 minutes** on the vibrator.)

Duplication

- Duplication can be carried out in gel moulds or in silicone moulds.
When working with a pressure compaction unit, silicone moulds and the duplicate model must be made under the same conditions (2 – 4 bar)! Duplicate in gel moulds only without pressure!
- Fill duplication mould on the vibrator and then remove it immediately from the vibrator.

Mixing ratio 100 g Powder : 20 ml Mixing liquid

for 2 duplicate models


	WiroFine	Liquid	Aqua dest.	Mixing liquid	
				Total	Concentration
Liquid: BegoSol K	1 x 400 g	56 ml	24 ml	80 ml	70%

	Gel (Castogel®, Wirodouble®, WiroGel M)	Silicone (Wirosil®)
Remove master model	after 45 – 60 minutes	after 45 – 60 minutes
Surface treatment	<i>Durol</i>	<i>Durofluid*</i>
Drying	30 minutes (250 °C / 500 °F)	10 minutes (80 – 100 °C / 180 – 210 °F)
Dipping/spraying	briefly 3 times (approx. 2 seconds)	Do not exceed 100°C/210°F
Subsequent drying	5 minutes (250 °C / 500 °F)	5 minutes (80 – 100 °C / 180 – 210 °F)

* do not use for duplication with gel

After removal of master model the silicone mould MUST bench set for 20-30 minute before mfg of refractory model!

Investment

- Before investing, prepare the wax-up with *Wiropaint plus* fine investment material or *Aurofilm* wetting agent (please follow the processing instructions for the products).
- Fill mould ring on the vibrator and then take away from vibrator immediately.
- Remove the mould ring 10 minutes after investing with a total bench set of 20 minutes (even if under pressure).

Mixing


Mixing ratio 100 g Powder : 20 ml Mixing liquid

for 1 mould

	WiroFine	Liquid	Aqua dest.	Mixing liquid	
				Total	Concentration
Liquid: BegoSol K	1 x 400 g	56 ml	24 ml	80 ml	70%

Precision casting investment material for partial dentures – shock or conventional heating

Preheating

	Conventional preheating	Shock heat
Setting after investment	at least 30 minutes	20 minutes
Insertion temperature	room temperature	700 – 1000 °C / 1290 – 1830 °F
Holding levels	250 °C / 500 °F (5 °C/min / 9 °F/min) 570 °C / 1060 °F (7 °C/min / 12 °F/min)	– –
Final temperature	950 °C – 1050 °C / 1740 °F – 1920 °F	
Holding times	30 – 60 minutes (depending on size and number of moulds)	


Place sprue former face down while preheating!

Risk of injury in connection with shock heating. Caution: Danger of darting flame!

Place all moulds in the furnace within 10 seconds – then keep the furnace door closed for 15 minutes!

After casting

After casting allow the moulds to cool down until warm to the touch, **do not quench in water!**

Investment materials contains quartz. Do not inhale dust! Danger of lungs harms (silicosis, lung cancer).

To avoid dust during deflasking, place the moulds in water after they have cooled down completely after casting until they are thoroughly moistened.

Data

Processing time at 20 °C / 70 °F approx. 3.5 min

Characteristic material values in accordance with DIN EN ISO 9694

(70 % BegoSol® K)

Beginning of setting (Vicat time) approx. 6.0 min

Compressive strength


(after 2 hours) approx. 11 MPa

Linear thermal expansion approx. 0.8 %


Flowability approx. 130 mm

This product was made according to the specifications of DIN EN ISO 9694 and meets its requirements in all respects.

Thermal expansion curve WiroFine

Availability and recommendations

WiroFine	1 carton 18 kg = 45 400 g bags	– Order No. 54345
	1 carton 6 kg = 15 400 g bags	– Order No. 54344
BegoSol® K	1 bottle = 1000 ml	– Order No. 51120
	1 canister = 5000 ml	– Order No. 51121
BegoSol®	1 bottle = 1000 ml	– Order No. 51090
	1 canister = 5000 ml	– Order No. 51091

<i>EasyMix</i>	26090 (230 V)	<i>Wirosil®</i>	52001 (2 kg)	<i>DuroI</i>	52111 (1000 ml)
<i>Castogel®</i>	52052 (6 kg)	<i>Wirosil® duplicating flask system</i>		<i>Durofluid</i>	52008 (100 ml)
<i>Wirodouble®</i>	52050 (6 kg)		52072 (small)	<i>Wiropaint plus</i>	51100 (200 ml)
<i>WiroGel M</i>	54351 (6 kg)		52083 (large)	<i>Aurofilm</i>	52019 (100 ml)

Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products are subject to a constant further development. Therefore alterations in construction and composition are reserved.

For particularly good results we recommend an alloy from the following groups, depending on the indication

Tel. +49 421 2028-282
www.bego.com

WIRONIUM®
(Only given to authorized laboratories)

Wironit®

PlatinLloyd®

since 1890